

Lab Report of Laser Spectroscopy Demonstration of 200kHz System

Tao LI

September 30, 2015

1 Questions

1. What do we mean by intracavity frequency doubling in the pump laser?

It means that a non-linear crystal, which could generate doubling frequency under SHG, is placed within the laser cavity. And a dichroic mirror in the cavity could only permit doubling frequency out while has a high reflection index to fundamental frequency.

2. What is the mode locking mechanism used in the Ti:sapphire laser?

Ti:sapphire crystal has Kerr effect, a non-linear effect of refractive index that $n_m = n_0 + I n_2$. When placed in the resonator, it has no(slight) focusing effect of weak light while having strong focusing effect of intense light. After travelling back and forth for many times, ultimately it would produce only a single, intense pulse. This is called *passive mode locking*.

3. Why do we need two prisms in the cavity?

The two prisms form a prism compressor system that to compensate the intra-cavity laser dispersion of Ti:sapphire laser. In detail, when laser beam goes through the first prism, different wavelengths would travel different paths. After the second prism, all wavelengths are adjusted to the same direction with different time, which would shorten the duration of laser pulse.

4. What is the function of the adjustable slit in the Ti:sapphire laser?

Adjustment of slit width changes the bandwidth of laser and inhibits secondary laser peak.

5. What is a regenerative amplifier?

In the regenerative amplifier, the input seed pulse enters the amplifier cavity and the Ti:sapphire crystal and a short and intense pulse is generated. After multiple passes in the gain crystal, the amplified pulse exits the amplifier.

6. What type of Q-switch is used in the amplifier? What is the function of the Q-switcher?

Electro-optic Q-switch.

It requires a linear polarization filter and an EO polarization rotator in the resonator. After low intense laser being rotated of polarization, the beam cannot pass through the linear polarization filter. Until it develops to be a giant and intense laser, the beam outputs when the linear polarization filter is removed. In short, EO Q-switch is to generate high intense laser beam.

7. What method is used to compress the pulse in the system?

Mode locking with Kerr effect is used to compress the laser to short pulse and a pair of prisms is used to shorten the pulse duration by compensation of temporal dispersion.

Consequently, chirped pulse amplification (CPA) also compresses the dispersed pulse.

8. Why do we need to compress the pulses?

We compress the pulse to a short duration pulse that helps to study fast chemical reaction and generate extremely high power.

9. What is the range of wavelength available for the OPA(Signal and Idler)? What is the wavelength of the idler if the signal is 580nm?

The wavelength available for signal and idler should obey the momentum conservation principle and phase matching condition. In general, $n_p v_p = n_i v_i + n_s v_s$.

Suppose the pump beam has wavelength of 500nm, according to

$$\frac{1}{\lambda_p} = \frac{1}{\lambda_i} + \frac{1}{\lambda_s}$$

we find the wavelength of idler is $4.67\mu m$.